

Negative pressure transport chamber for infected patients



Airtight chamber allowing infected patients to be transported under negative pressure.

The device includes:

- Chamber made from 300 micron crystal PU
- Airtight zip fastener allowing chamber to be opened and closed
- Three P3 inlet filters allowing filtration of outside air and preventing air from being released in case of occasional overpressure
- Two P3 air intake filters linked to the vacuum motor
- Vacuum motor with 160l/min capacity and 8h autonomy
- 2 gloveboxes allowing handling of the patient
- 2 airtight plugs for tubes (drips/high concentration mask)
- 3 patient straps
- 2 stretcher straps
- Carry bag

Dimensions of chamber:

- **Length:** 2060mm
- **Width:** 750mm
- **Height:** 510mm
- **Weight:** 7kg

Optional extras and accessories on request



الرواد

رواد للخدمات الطبية ROWAD MEDICAL SERVICES



4710 Prince Mohammed Bin Abdulaziz RD
Al Sulaimaniyah Dist - Riyadh 12243 - 6934
Saudi Arabia



Info@alrowad-ksa.com



Tel.: +966 114 005 640
Mob.: +966 582 748 097

OUR BRANCHES



SAUDI ARABIA



UNITED ARAB EMIRATES



SWITZERLAND