



 surewash

YOUR PARTNER IN HAND HYGIENE TRAINING

SureWash
GO



SYSTEM BENEFITS

SureWash is a cost-effective hand hygiene training system teaching the WHO technique.



200+

hospitals groups and universities use the SureWash platform




650,000+

Registered users on SureWash to date



1,000,000+

hand hygiene training sessions are delivered by SureWash annually



SUREWASH GO

The SureWash GO is a portable training system that uses a live video camera to measure hand motions. It provides real-time feedback to support the user in a “Deliberate Practice” learning framework.

This approach helps users to develop muscle memory, therefore improving compliance, reducing the risk of infections and improving patient outcomes.

Not only is this training system beneficial to hospital staff, SureWash GO also generates awareness around infection control amongst patients and family members. Its ease of use results in patients and visitors finding the training an engaging and positive experience.

As it weighs just 4.5kg, the GO device is fully portable and comes with its own molded case that has custom foam for ease of transport.

Ultimately, SureWash not only saves organisations time, money and productivity, but has a very real impact in improving patient outcomes.



How effective is SureWash?



Improving long term results on hand hygiene interventions


REPORTING SUITE

The Reporting Suite on SureWash.NET is an audited training record for regulatory bodies.

The reports allow for the automation and aggregation of data across organisations. Data can be defined by role or department and can be narrowed to a certain date. Management reports can also be downloaded and printed in PDF or Excel format.

SureWash.NET has a quiz option which provides instant feedback regarding the standard of knowledge in the organisation while also helping to address the gaps in learning. Along with this, lessons are available on PowerPoint presentation to deliver sharp bursts of learning. This content can be personalised by the organisation to encourage engagement and increased interaction among users.

All elements of the system combine to provide a complete solution for hand hygiene training.




WHAT OUR USERS SAY

“Not another boring competency. SureWash uses an interactive video game system to improve hand hygiene skills. Staff love competing with each other and don’t even realize they are learning proper hand hygiene technique.”

Nancy Vallande
Miriam Hospital, Providence

“I no longer need to think about the order of the WHO 6 pose technique because it's in my memory from repeating it multiple times on the SureWash machine.”

Lenora Leonard
The Beacon Hospital, Dublin



Specifications

Height: From 43cm/ 17 inches for wall mounting to 64cm/25 inches with camera up in tabletop mode.

Base Dimensions: 4.7 x 10.7 inches / 12 x 27cm

Weight: 4.5 k Kg / 9.92 lb (10 kg / 22 lb with Case)

Portable Carry Case: Pele Storm 5100 with custom foam

Components: Microsoft Surface Pro (with Intel Core i5).
SoftKinetic DS325 3D Camera.
WiFi.
External ports: 1 x USB 2.0.
Power 110 / 240 volt.
Kensington security lock slot.

RFID Card Reader.



رواد للخدمات الطبية
ROWAD MEDICAL SERVICES

www.alrowad-ksa.com
Al sulaimaniyah Dist - Riyadh
Saudi Arabia
info@alrowad-ksa.com
mob: +966 58 274 8097